

Cell and Tissue Culture Flasks , Non-treated

| Cat. No.             | Volume (mL) | Surface Type | Recommended Working Volume (mL) | Type of Cap | Sterile | Qty. Per Bag/Case |
|----------------------|-------------|--------------|---------------------------------|-------------|---------|-------------------|
| TCF001025            | 25          | Non-treated  | 8                               | Plug Seal   | Y       | 10/200            |
| TCF002025            | 25          | Non-treated  | 8                               | Vent        | Y       | 10/200            |
| TCF001050            | 50          | Non-treated  | 17.5                            | Plug Seal   | Y       | 10/200            |
| TCF002050            | 50          | Non-treated  | 17.5                            | Vent        | Y       | 10/200            |
| TCF001250            | 250         | Non-treated  | 60                              | Plug Seal   | Y       | 5/100             |
| TCF002250            | 250         | Non-treated  | 60                              | Vent        | Y       | 5/100             |
| TCF001600            | 600         | Non-treated  | 125                             | Plug Seal   | Y       | 5/40              |
| TCF002600            | 600         | Non-treated  | 125                             | Vent        | Y       | 5/40              |
| TCF101600 (Extended) | 600         | Non-treated  | 200                             | Plug Seal   | Y       | 5/40              |
| TCF102600 (Extended) | 600         | Non-treated  | 200                             | Vent        | Y       | 5/40              |
| TCF001225            | 850         | Non-treated  | 200                             | Plug Seal   | Y       | 5/25              |
| TCF002225            | 850         | Non-treated  | 200                             | Vent        | Y       | 5/25              |
| TCF001850            | 850         | Non-treated  | 200                             | Plug Seal   | Y       | 3/18              |
| TCF002850            | 850         | Non-treated  | 200                             | Vent        | Y       | 3/18              |

Cell and Tissue Culture Flasks, TC treated

| Cat. No.             | Volume (mL) | Surface Type | Cell Culture Surface Area (cm²) | Recommended Working Volume (mL) | Type of Cap | Sterile | Qty. Per Bag/Case |
|----------------------|-------------|--------------|---------------------------------|---------------------------------|-------------|---------|-------------------|
| TCF011025            | 25          | TC treated   | 12.5                            | 8                               | Plug Seal   | Y       | 10/200            |
| TCF012025            | 25          | TC treated   | 12.5                            | 8                               | Vent        | Y       | 10/200            |
| TCF011050            | 50          | TC treated   | 25                              | 17.5                            | Plug Seal   | Y       | 10/200            |
| TCF012050            | 50          | TC treated   | 25                              | 17.5                            | Vent        | Y       | 10/200            |
| TCF011250            | 250         | TC treated   | 75                              | 60                              | Plug Seal   | Y       | 5/100             |
| TCF012250            | 250         | TC treated   | 75                              | 60                              | Vent        | Y       | 5/100             |
| TCF011600            | 600         | TC treated   | 182                             | 125                             | Plug Seal   | Y       | 5/40              |
| TCF012600            | 600         | TC treated   | 182                             | 125                             | Vent        | Y       | 5/40              |
| TCF111600 (Extended) | 600         | TC treated   | 182                             | 200                             | Plug Seal   | Y       | 5/40              |
| TCF112600 (Extended) | 600         | TC treated   | 182                             | 200                             | Vent Cap    | Y       | 5/40              |

| Cat. No.  | Volume (mL) | Surface Type | Recommended Working Volume (mL) | Type of Cap | Sterile   | Qty. Per Bag/Case |      |
|-----------|-------------|--------------|---------------------------------|-------------|-----------|-------------------|------|
| TCF011225 | 850         | TC treated   | 225                             | 200         | Plug Seal | Y                 | 5/25 |
| TCF012225 | 850         | TC treated   | 225                             | 200         | Vent      | Y                 | 5/25 |
| TCF011850 | 850         | TC treated   | 300                             | 200         | Plug Seal | Y                 | 3/18 |
| TCF012850 | 850         | TC treated   | 300                             | 200         | Vent      | Y                 | 3/18 |

Cell and Tissue Culture Flasks, CellATTACH<sup>®</sup> treated

| Cat. No.            | Volume (mL) | Surface Type                    | Cell Culture Surface Area (cm <sup>2</sup> ) | Recommended Working Volume (mL) | Cap Style | Sterile | Qty. Per Bag/Case |
|---------------------|-------------|---------------------------------|----------------------------------------------|---------------------------------|-----------|---------|-------------------|
| CAF011025           | 25          | CellATTACH <sup>®</sup> treated | 12.5                                         | 8                               | Plug Seal | Y       | 10/200            |
| CAF012025           | 25          | CellATTACH <sup>®</sup> treated | 12.5                                         | 8                               | Vent      | Y       | 10/200            |
| CAF011050           | 50          | CellATTACH <sup>®</sup> treated | 25                                           | 17.5                            | Plug Seal | Y       | 10/200            |
| CAF012050           | 50          | CellATTACH <sup>®</sup> treated | 25                                           | 17.5                            | Vent      | Y       | 10/200            |
| CAF011250           | 250         | CellATTACH <sup>®</sup> treated | 75                                           | 60                              | Plug Seal | Y       | 5/100             |
| CAF012250           | 250         | CellATTACH <sup>®</sup> treated | 75                                           | 60                              | Vent      | Y       | 5/100             |
| CAF011600           | 600         | CellATTACH <sup>®</sup> treated | 182                                          | 200                             | Plug Seal | Y       | 5/40              |
| CAF012600           | 600         | CellATTACH <sup>®</sup> treated | 182                                          | 200                             | Vent      | Y       | 5/40              |
| CAF111600(Extended) | 600         | CellATTACH <sup>®</sup> treated | 182                                          | 200                             | Plug Seal | Y       | 5/40              |
| CAF112600(Extended) | 600         | CellATTACH <sup>®</sup> treated | 182                                          | 200                             | Vent      | Y       | 5/40              |
| CAF011850           | 850         | CellATTACH <sup>®</sup> treated | 300                                          | 200                             | Plug Seal | Y       | 3/18              |
| CAF012850           | 850         | CellATTACH <sup>®</sup> treated | 300                                          | 200                             | Vent      | Y       | 3/18              |