

Syringe Filters

Cat. No.	Membrane Material	Color	Pore Size(μm)	Housing Diameter (mm)	Sterile	Qty. Per Box/ Case
FMC201013	MCE	Yellow	0.22	13	Y	100/800
FMC201025			0.22	25	Y	45/360
FMC201030			0.22	30	Y	45/360
FMC401013			0.45	13	Y	100/800
FMC401025			0.45	25	Y	45/360
FMC401030			0.45	30	Y	45/360
FPV103013	PVDF	Blue	0.1	13	Y	100/800
FPV103025			0.1	25	Y	45/360
FPV103030			0.1	30	Y	45/360
FPV203013			0.22	13	Y	100/800
FPV203025			0.22	25	Y	45/360
FPV203030			0.22	30	Y	45/360
FPV403013			0.45	13	Y	100/800
FPV403025			0.45	25	Y	45/360
FPV403030			0.45	30	Y	45/360
PTF205013	PTFE	White	0.22	13	Y	100/800
PTF205030			0.22	25	Y	45/360
PTF205030			0.22	30	Y	45/360
PTF405013			0.45	13	Y	100/800
PTF405025			0.45	25	Y	45/360
PTF405030			0.45	30	Y	45/360
FNY202013	Nylon	Pink	0.22	13	Y	100/800
FNY202025			0.22	25	Y	45/360
FNY202030			0.22	30	Y	45/360
FNY402013			0.45	13	Y	100/800
FNY402025			0.45	25	Y	45/360
FNY402030			0.45	30	Y	45/360
FPE204013	PES	Green	0.22	13	Y	100/800
FPE204025			0.22	25	Y	45/360
FPE204030			0.22	30	Y	45/360
FPE404013			0.45	13	Y	100/800
FPE404025			0.45	25	Y	45/360
FPE404030			0.45	30	Y	45/360
FCA206013	CA	Red	0.22	13	Y	100/800
FCA206025			0.22	25	Y	45/360
FCA206030			0.22	30	Y	45/360
FCA406013			0.45	13	Y	100/800
FCA406025			0.45	25	Y	45/360
FCA406030			0.45	30	Y	45/360
SCA207013	SFCA	Deep Bule	0.22	13	Y	100/800
SCA207025			0.22	25	Y	45/360
SCA207030			0.22	30	Y	45/360
SCA407013			0.45	13	Y	100/800
SCA407025			0.45	25	Y	45/360
SCA407030			0.45	30	Y	45/360
FPE204113			0.22	13	Y	100/800
FPE204125			0.22	25	Y	45/360

Cat. No.	Membrane Material	Color	Pore Size(μm)	Housing Diameter (mm)	Sterile	Qty. Per Box/ Case
FPE204130	PES Express	Green	0.22	30	Y	45/360
FPE404113			0.45	13	Y	100/800
FPE404125			0.45	25	Y	45/360
FPE404130			0.45	30	Y	45/360
GFA201025	GF1.1um+C A0.22um	Natural	0.22	25	Y	45/360
GFA201030			0.22	30	Y	45/360
GFA401025			0.45	25	Y	45/360
GFA401030			0.45	30	Y	45/360