

MS-S

Magnetic Stirrer



Specifications	MS-S
Work plate Dimension [W x D]	Ø135mm (5 inch)
Work plate material	Stainless steel cover with ceramic
Motor type	Brushless DC motor
Stirring positions	1
Max. stirring quantity[H ₂ O]	20L
Max.magnetic bar[length]	80mm
Speed range	0-1500rpm
Speed display	Scale
Protection class	IP21
Motor input Power	18W
Motor output Power	10W
Power	30W
Voltage,Frequency	100-240V, 50/60Hz
Dimension[WxDxH]	160x280x85mm
Weight	2.8kg
Permissible ambient temperature and humidity	5-40°C, 80%RH

